

# FORMFIT® WALKER

## INDICATIONS FOR USE

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative, trauma or rehab use

## FEATURES AND BENEFITS

- High and low top versions available
- Shock absorbing sole reduces impact during heel strike
- Low, wide rocker bottom promotes a stable, natural gait
- Breathable, seamless liner

## REQUIRED MEASUREMENTS

- Shoe size



## ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Standard	Low Top
X-Small	2 – 4	3.5 – 5.5	1.5 – 3.5	33 – 35.5	W0200BLK	W0300BLK
Small	4.5 – 7	6 – 8	4 – 6.5	36 – 39	W0400BLK	W0500BLK
Medium	7.5 – 10.5	8.5 – 11.5	7 – 10	39.5 – 44	W0600BLK	W0700BLK
Large	10.5 – 12.5	11.5 – 13.5	10.5 – 12	44.5 – 46.5	W0800BLK	W0900BLK
X-Large	12.5 +	13.5 +	12.5 – 14.5	47 – 50	W1000BLK	W1100BLK

## ACCESSORIES

### REPLACEMENT SOFTGOODS

Size	Standard	Low Top
X-Small	WKSG200BLK	WKSG0300BLK
Small	WKSG400BLK	WKSG500BLK
Medium	WKSG600BLK	WKSG700BLK
Large	WKSG800BLK	WKSG900BLK
X-Large	WKSG1000BLK	WKSG1100BLK

### OTHER ACCESSORIES

Description	Part#
10 cm softgood width extension	WKSG-1BLK
Tibia and calf plates for use when additional splinting is required	WKSP100
Hot/Cold therapy insert for safe and reusable cryotherapy	WKGEL100
Heel wedge, Small	AWK0001
Heel wedge, Medium	AWK0005
Heel wedge, Large	AWK0009
Adhesive insole for heel wedge	WHW-W0107
Walker Strap Extender	10882